

SURGICAL BLADES




Small Fitment Ho.3



Large Fitment Ho.4

SURGICAL BLADES

Available in both stainless and carbon steel. Numbers 10, 10A, 11, 11P, 12, 12B/12D, 13, 14, 15, 15T, 15C, 15A, 16, 17, 40 fit handle numbers 3, 3L, 5, 7 & 9. Numbers 18, 19, 20, 21, 22A, 22, 23, 24, 25, 26, 36 fit handle numbers 4, 4L & 6. Individually V.C.I paper and foil packed. Supplied sterile and non sterile, in box packing of 100.



Surgical Handle No.3



Surgical Handle No.4



DISPOSABLE SCALPELS

Available in pattern numbers

10, 10A, 11, 11P, 12, 12B, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26 supplied sterile. Stainless Steel blade tempered with plastic handle. Protected by polythene blade guard. Individually packed in peel pack medical grade paper. 10 pieces in a box. Also available non sterile bulk pack 100 pcs in a box.



THUMB MINI SCALPELS

Can be used as a part of the "Ward Kit". Small plastic handles, supplied non-sterile, fitted with blade numbers 10, 11, 15, 22 & 60. Protected with a polythene blade guard. Bulk packed 100 pieces in a box.

