



EN

Extracting Forceps

DE

Extraktionszange

ES



Fórceps de extracción

Page # 01 To 28



Handle Selection For Extracting Forceps

Dear Valued Buyers if any Universal Fig. No. of Extracting Forceps not Found in our Catalogue you Can Refer the Fig. No. From These Diagram.



Dear Valued Buyers these are our Handle Styles



English Pattern



English Pattern

English Pattern,
Anatomical Handle

American Pattern

Deep Grip Forceps
With Hole Pattern

Hole Pattern

**A****B****C****D****E****F**

Function: To remove teeth from bony socket



DE-001
Fig. 1
Upper Centrals
and Canines



DE-002
Fig. 2
Upper Laterals
and Canines



DE-003
Fig. 4
Lower Incisors
and Canines



DE-004
Fig. 7
Upper Premolars



DE-005
Fig. 8
Lower Premolars



DE-006
Fig. 13
Lower Premolars



DE-007
Fig. 17
Upper Molars,
Right



DE-008
Fig. 18
Upper Molars, Left



DE-009
Fig. 21
Lower Molar



DE-010
Fig. 22
Lower Molar



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Function: To remove teeth from bony socket



DE-011
18A
upper molars



DE-012
19
upper third molars



DE-013
20
lower third molar



DE-014
24
lower molar, right



DE-015
29
upper roots



DE-016
30
Upper Roots



DE-017
33
lower roots



DE-018
33A
Lower Roots



DE-019
33M
Lower roots
with long beak



DE-020
46L



DE-021
51
Upper Roots



DE-022
51A
Upper Roots



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines



Function: To remove teeth from bony socket

WISDOM EXTRACTING FORCEPS



DE-023
Routerier
(# 22 1/2 L)
lower molars and
wisdoms, left



DE-024
Routerier
(# 22 1/2 R)
lower molars and
wisdoms, right



DE-025
Routerier
(# 33 1/2 L)
lower molar and
wisdom roots, left



DE-026
Routerier
(# 33 1/2 R)
lower molar and
wisdom roots, right



*Routerier extracting forceps for lower molars
and third molars, right*



DE-027
Fig. 186L
lower molars



DE-028
Fig. 186R
lower molars



DE-029
51AL



DE-030
67
Upper Third Molars



DE-031
67A
Upper Third Molars



DE-032
79
Lower Third Molars





DE-033
73
Lower Molars



DE-034
73 A
lower molars



DE-035
73 S



DE-036
74
Lower Roots



DE-037
74N
Lower Roots



DE-038
75
Lower Premolars



DE-039
86A
Lower Molars



DE-040
86B
lower molars



DE-041
86C
lower molars



DE-042
87
Lower Molars



DE-043
89
Upper Molars,
Right



DE-044
90
Upper Molars,
Left



Straight extracting forceps for upper incisors and canines




Straight extracting forceps for upper incisors and canines



Straight extracting forceps for upper incisors and canines

A-TRAUMATIC EXTRACTION
 improved tapered beaks for complete root contact. Cross serrations. Deep Gripe Extracting Forceps



• No Fracture • No Tissue Damage • Extraction With Easy



Thin Tip

DE-045
34

Upper incisors and
 premolars
 deep-gripping



Wide Tip

DE-046
34A

Upper centrals / cains
 deep-gripping



Thin Tip

DE-047
Fig. 35

upper permolars
 deep-gripping



Wide Tip

DE-048
35A

Upper incisors and
 premolars
 deep-gripping



Thin Tip

DE-049
36

Lower premolars
 deep-gripping



Wide Tip

DE-050
36A

Lower premolars
 deep-gripping



DE-050
336

lower permolars
 deep-gripping



DE-052
F-7C

lower molar



DE-053
9

Apex Forceps Lower Universal
 Deep Root Pick



DE-054
79N

Lower third molars
 deep-gripping



DE-055
79

Lower Thrid Molars



Function:
 Extracting forceps are used to further loosen teeth by rotational and tilting movements and finally to extract the teeth in combination with a fulcrum.



DEEP ROOT PICK FORCEPS



DE-056
146

Apex Forceps Lower Root
Deep Root Pick



DE-057
153

Apex Forceps Upper Root
Deep Root Pick



DE-058
149

Apex Forceps Upper Root
Deep Root Pick



DE-061
Fig. 44



DE-060
Fig. 44Z



DE-059
Fig. 45Z



DE-062
46XL
Apex Forceps Lower Root
Deep Root Pick



DE-063
67
Upper Third Molars



DE-064
67A
Upper Third Molars



forceps for upper incisors, premolars and roots



DE-065
67N
Lower third molars
deep-gripping



DE-066
331
upper roots
deep-gripping



Fig. 1
upper
laterals
and canines



DE-067
Set of 10 Pieces
Fig. 1
Fig. 7
Fig. 13
Fig. 17
Fig. 18
Fig. 22
Fig. 51A
Fig. 67
Fig. 74N
Fig. 79

Fig. 7
upper premolars



Fig. 13
lower premolars



Fig. 17
upper molars, right



Fig. 18
upper molars, left



Fig. 22
lower roots



Fig. 51A
upper roots



Fig. 67
upper third molars



Fig. 74N
lower molars



Fig. 79
lower third molars



Extracting forceps for upper roots





Pediatric Tooth Extracting Forceps Set



DE-069
English Pattern Klein #51S
Pedodontic Upper root



DE-068
English Pattern Klein #137
Pedodontic, Upper Incisors Spring handle



DE-072
English Pattern Klein #5
Pedodontic Lower Incisors



DE-070
English Pattern Klein #139
Pedodontic Upper Premolars

Baby Tooth Extraction Forceps




DE-075
Fig. # 01 to 07
Set of 7 Pieces



DE-073
English Pattern Klein #7
Pedodontic Lower Roots



DE-071
English Pattern Klein #3
Pedodontic Upper Molars Universal



DE-074
English Pattern Klein #6
Pedodontic Lower Molars, Spring Handle



**Witzel
DE-076**
170mm
Upper roots, Universal



**Witzel
DE-077**
165mm
Lower roots, Universal



DE-078
#400(300)
165mm
Upper roots, Universal



DE-079
#401(301)
160mm
Lower roots, Universal



DE-080
Witzel Forceps With Gold
Set of 3 Pcs



DE-084
Witzel Forceps With Hole And Screw Pattern
Set of 3 Pcs



DE-088
Witzel Forceps With Anatomical Handle
Set of 3 Pcs



DE-081
#401(301)
160mm
Lower roots, Universal



DE-085
#401(301)
160mm
Lower roots, Universal



DE-089
#401(301)
160mm
Lower roots, Universal



DE-082
#400(300)
165mm
Upper roots, Universal



DE-086
#400(300)
165mm
Upper roots, Universal



DE-090
#400(300)
165mm
Upper roots, Universal



**Witzel
DE-083**
165mm
Lower roots, Universal



**Witzel
DE-087**
165mm
Lower roots, Universal



**Witzel
DE-091**
165mm
Lower roots, Universal



DE-092
 Splinter Forceps
 15.5cm



Stieglitz Root Splinter Forceps



DE-093
 Splinter Forceps
 15.5cm 45°



DE-094
 Stieglitz
 Root Fragments




DE-095
 Stieglitz
 Root Fragments



DE-096
 Stieglitz Splinter Forceps
 140mm 90°



DE-097
 Stieglitz Splinter Forceps
 140mm Straight




With Diamond

DE-098
130mm



With Diamond

DE-099
130mm



With Diamond


DE-100
125mm



DE-101
130mm



DE-102
130mm



DE-103
125mm



DE-104
#1STD
Winter
upper centrals and canines



DE-105
Fig. 1
upper incisors
and roots



DE-106
Cryer
Fig. 150



DE-108
 Fig. 62
 incisors, premolars,
 deciduous teeth



DE-109
 # 10S
 upper roots



DE-110
 # 17
 lower molars



DE-111
 Harris
 Fig. 18R



DE-112
 Harris
 Fig. 18L



DE-107
 Fig. 150A
 Cryer
 upper incisors, premolars, roots



DE-113
Fig. 210
upper third molars



DE-114
Fig. 53R
upper molars, right



DE-115
Fig. 53L
upper molars, left



DE-116
NEVIUS
Fig. 88R





DE-117
Fig. 24
upper molars



DE-118
210S
upper third molars



DE-119
210H
upper third molars



DE-120
286
Upper Roots,
Premolars, Incisors



DE-121
99C
Kells
upper incisors,
canines, premolars




DE-122
222
lower third molars




DE-123
286
Lower Molar




DE-124
Nevius
Fig. 88L
upper molars, left



DE-125
Fig. 65
upper incisors
and roots



DE-126
Parmlly Fig. 32A
upper premolars
and roots



DE-127
Parmlly Fig. 32
upper premolars
and molars



DE-128
 COHEN
 universal

DE-129
 Fig. 203
 lower incisors
 and roots

DE-130
 Tomes Fig. 69
 upper and lower
 roots



DE-131
Fig. 23
lower molars

DE-132
Fig. 16
lower molars



DE-133
Physic Fig. 5
lower third molars



DE-134
Harris Fig. 15
lower molars



DE-135
Fig. 103




DE-136
Hull
Fig. 101
 Deciduous teeth



DE-137
Cryer
Fig. 150S
 upper teeth,
 for children




DE-138
Fig. 150 AS



DE-139
Cryer
 Fig. 151
 lower teeth,
 for children

DE-140
 Fig. 151 S
 lower incisors,
 premolars, roots


DE-141
 Fig. 151 AS




DE-142
16S
 lower molars for children



DE-143
17S
 lower molars



DE-144
23 S
 Lower molars




DE-145
27
 lower molars




DE-146
18LS
 upper molars left



DE-147
18RS
 upper molars right




DE-148
85 Tomas
 lower canines, premolars and molars




DE-149
85A Tomas
 lower canines, premolars and molars



DE-150
69S
 upper and lower roots



DE-151
99A Kells
 upper incisors, canines, premolars



DE-152
217
 lower molars



DE-153
#10 SK
Upper Molars Sup



DE-154
#17 SK Universal
Lower Molars



DE-155
#23 SK Universal
Upper Molars



DE-156
#69 SK Tomas
upper and lower roots



DE-157
#99 SK
upper incisors,
canines, premolars



DE-158
#150 SK
upper incisors,
premolars, roots



DE-159
#151 SK
lower incisors,
premolars, roots



DE-160
#222 SK
lower third molars



DE-161
Fig. 1
Upper Anterior

DE-162
Lower Molar
Fig. 23 Cow Horn



DE-163
Fig. 74
Lower Anterior



DE-164
Fig. 88L
Upper Molar Left



DE-165
Fig. 88R
Upper Molar Right



DE-166
Fig. 74
Lower Anterior



DE-167
Fig. 106
Lower incisors



DE-168
Fig. 174
Extra Narrow
Lower Anterior
Apical Retention
Forceps



DE-169
Fig. 859
Lower Molars
Roots





Carbide Diamond Coated Point



DE-170
150 T.C.
Carbide Serrated Beak
Upper Incisors, Cuspids,
Bicusps & Roots,
Universal



DE-171
151 T.C.
Carbide Serrated Beak
Lower Incisors, Cuspids,
Bicusps & Roots, Universal



DE-172
1STD T.C.
Carbide Serrated Beak
Upper Incisors, Cuspids,



DE-173
17 T.C.
Carbide Serrated Beak
1st & 2nd Lower Molars,
Universal



DE-174
53L T.C.
Carbide Serrated Beak
1st & 2nd Upper Molars,
Left



DE-175
53R T.C.
Carbide Serrated Beak
1st & 2nd Upper Molars,
Right



DE-176
65 T.C.
Carbide Serrated Beak
Upper Incisors ,Roots &
Fragments



DE-177
69 T.C.
Carbide Serrated Beak
Upper & Lower Fragments
And Roots



DE-178
150S T.C. Pedro forcep
Carbide Serrated Beak
Pedo, Upper Primary Teeth
and Roots, Universal



DE-179
151S T.C. Pedro forcep
Carbide Serrated Beak
Pedo, Lower Primary Teeth
and Roots, Universal





DE-180
Blumenthal
15cm



DE-181
Blumenthal
15cm



DE-182
Blumenthal
15cm



DE-183
Friedman
14cm



DE-184
Friedman
14.5cm



DE-185
Micro-Friedman
14cm



DE-186
Micro-Friedman
15cm



DE-187
Micro-Friedman
16cm



DE-188
Kleinert-Kutz
18.5cm



DE-189
Beyer
18.5cm



DE-190
Mead 1A
18.5cm



DE-191
Luer
18.5cm



DE-192
Meissner
Extracting Forceps

DE-193
Meissner
Extracting Forceps

DE-194
Meissner
Extracting Forceps

DE-195
Comellas Forceps
Anatomical Handle



DE-196 1.8ml
DE-197 2.2ml



DE-198 1.8ml
DE-199 2.2ml



DE-200 1.8ml
DE-201 2.2ml



DE-202 1.8ml
DE-203 2.2ml



DE-204 1.8ml
DE-205 2.2ml



DE-206 1.8ml
DE-207 2.2ml




DE-208
1001
Intraligamental
Anesthesia
Syringe in
Black Colour



DE-209 1.8ml
DE-210 2.2ml
Angled




DE-211
Aluminium - Dark 1.8 ml



DE-212 Aluminium - Silver 1.8 ml
DE-213 Aluminium - Orange 1.8 ml
DE-214 Aluminium - Green 1.8 ml
DE-215 Aluminium - Blue 1.8 ml



Type C
DE-216 1,8ml
DE-217 2,2ml



Type C
DE-218 1,8ml
DE-219 2,2ml



2-Ring
DE-220 1.8ml
DE-221 2.2ml



Aspirating Syringe
DE-222 1.8mm
DE-223 2.2mm



Catridge Syringe
DE-224 1.8mm
DE-225 2.2mm



Catridge Syringe
DE-226 1.8mm
DE-227 2.2mm



Catridge Syringe
DE-228 1.8mm
DE-229 2.2mm



Euro Breech loading
DE-230 1,8ml
DE-231 2,2ml



EuroBreech loading
DE-232 1,8ml
DE-233 2,2ml



Manual Aspirating Syringe
DE-234 1,8ml
DE-235 2,2ml



Quatro
DE-236 1,8ml
DE-237 2,2ml



Trio
DE-238 1,8ml
DE-239 2,2ml



Swedish Type
DE-240 1,8ml
DE-241 2,2ml



Ject Ultima
DE-242 1,8ml
DE-243 2,2ml
Injects 0,05 ml with every click Supplied with
angled and straight tubes



Ject angled
DE-244 1,8ml
DE-245 2,2ml



Ject straight
DE-246 1,8ml
DE-247 2,2ml



Disc System
DE-248 1,8ml
DE-249 2,2ml



Disc System
DE-250 1,8ml

**DE-067**

Set of 10 Pieces

Fig. 2

Fig. 7

Fig. 13

Fig. 17

Fig. 18

Fig. 22

Fig. 33

Fig. 51A

Fig. 67

Fig. 79

Fig.2
upper laterals and caninesFig. 7
upper premolarsFig. 13
lower premolarsFig. 17
upper molars rightFig. 18
upper molars leftFig. 22
lower rootsFig. 33
upper rootsFig. 51A
upper rootsFig. 67
upper third molarsFig. 74
lower rootsFig. 79
lower third molars

Extracting forceps for upper roots





DE-259
 Set of 10 Pieses
 Fig. 2
 Fig. 7
 Fig. 13
 Fig. 17
 Fig. 18
 Fig. 22
 Fig. 33
 Fig. 51A
 Fig. 67
 Fig. 79

Fig. 2
 upper laterals and canines



Fig. 7
 upper premolars



Fig. 13
 lower premolars



Fig. 17
 upper molars, right



Fig. 18
 upper molars, left



Fig. 22
 lower roots



Fig. 33
 lower roots



Fig. 51A
 upper roots



Fig. 67
 upper third molars



Fig. 79
 lower third molars



Extracting forceps for upper roots



**DE-251**

Extracting Forceps
Set of 10 Pieces

Extracting Forceps Fig. 2
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79

**DE-252**

Extracting Forceps
Set of 10 Pieces

Extracting Forceps Fig. 2
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79

**DE-253**

Extracting Forceps
Set of 13 Pieces

Extracting Forceps Fig. 2
Extracting Forceps Fig. 7
Extracting Forceps Fig. 13
Extracting Forceps Fig. 17
Extracting Forceps Fig. 18
Extracting Forceps Fig. 22
Extracting Forceps Fig. 33
Extracting Forceps Fig. 51A
Extracting Forceps Fig. 67A
Extracting Forceps Fig. 79
Extracting Forceps Fig. 86C
Extracting Forceps Fig. 89
Extracting Forceps Fig. 90

**DE-254**

Extracting Forceps Children Pattern
Set of 10 Pieces

Extracting Forceps Fig. 13S
Extracting Forceps Fig. 22S
Extracting Forceps Fig. 30S
Extracting Forceps Fig. 33S
Extracting Forceps Fig. 37
Extracting Forceps Fig. 38
Extracting Forceps Fig. 39
Extracting Forceps Fig. 39L
Extracting Forceps Fig. 39R
Extracting Forceps Fig. 51S

**DE-255**

Extracting Forceps English adult Adult
Set of 10 Pcs

**DE-256**

Extracting Forceps American Adult
Set of 10 Pcs

**DE-257**


Extracting Forceps American Adult
Set of 8 Pcs

**DE-258**

Extracting Forceps Pedro
Set of 7 Pcs

YOUR DISTRIBUTOR:

HygeiaMedical



Hygeia Medical Ltd. Reg. 616490 EUDAMED SRN IE-AR-000023941 TRN 3957302JH

www.hygeia.icu



IRELAND - GERMANY - MALAYSIA - INDONESIA - PHILPPINES - AUSTRALIA