

# HygeiaMedical™

MALAYSIA EDITION 2025



SCAN QR TO  
DOWNLOAD PDF

## Hygeia Medical

HYGEIA MEDICAL IS MORE THAN A SUPPLIER OF MEDICAL INSTRUMENTS



### Expertise

Hygeia Medical has some of the most experienced and reputable experts in the field among its associates, partners and suppliers

### Experience

Hygeia Medical is working with people with combined centuries of experience in the fields of manufacturing, maintenance and repair of surgical instruments in the ASEAN region.

### Reliability

Hygeia Medical is being relied on by medical professionals from around the ASEAN region as their primary quality supplier of choice.

### Commitment

Hygeia Medical is committed to becoming the regional standard for affordable quality instruments and to “work global but act local” by sharing know-how and expertise with our local partners.

## Quality Management Systems

HYGEIA MEDICAL LTD AND OUR SUPPLIERS ARE ISO AND CE CERTIFIED



### Hygeia Medical Ltd

Hygeia Medical is MDR registered with EUDAMED [SRN IE-AR-000023941](#) and further ISO 9001:2015 (IE90692A thru 24/11/2026) and also ISO 13485:2016 (IE90692H thru 24/11/2026) as well as CE certified by LMS Certification Ltd of Manchester, England.



### Hygeia Instruments (Malaysia) Sdn Bhd

Hygeia Instruments is GDPMD(MYG11206674 thru 15/11/2026) certified for import, distribution, storage and documentation of medical devices by CARE Certification International (M) Sdn Bhd of Selangor, Malaysia. Concepts for additional certifications are being implemented step by step

### Third Party Suppliers

As a rule all our third party suppliers are ISO certified and compliant, hold the required CE declarations and provide full documentation for all deliveries. All instruments are again individually inspected by Hygeia Quality Control and lodged prior to shipping.



## HYGEIA MEDICAL LIMITED, EIRE

Kyleteologue House, 9 Cherrywood  
Dundalk 21222, Co. Louth, Ireland

Hygeia Medical Limited Ireland is the principal supplier of HygeiaMedical branded instruments, as well as of a number of other reputable brands from Germany, Turkey, Korea, Taiwan, etc. to partners and associates in the ASEAN region.

Originally named Clearview Europe Limited, the company was purchased by its current owners in 2017 for the express purpose to distribute general and proprietary surgical instruments in the ASEAN countries and offer technology transfer and medical finance innovations to associates and partners across the region.

## HYGEIA MEDIZINTECHNIK & -VERTRIEB

PF2513 Hindenburgstrasse  
46395 Bocholt, NRW, Germany

In 2019 one of the partners also established a presence in Germany (**Hygeia Medizintechnik**) in Bocholt, near Düsseldorf, to serve as liaison to German suppliers and establish distribution links for devices, components, and consumables sourced from manufacturers in ASEAN countries.

A second German presence which specializes in distribution of disposables and supplies the veterinary market, followed in 2021 (**Hygeia Medizinvertrieb**) in Wiesbaden, near Frankfurt.

## HYGEIA MEDICAL (AUSTRALIA) PTY LIMITED

Hygeia Medical and Instruments recently set up this marketing office with local partners in Australia to become the distributor of our brands.

## HYGEIA MEDICAL INSTRUMENTS PVT LTD

Hygeia Medical and Instruments established this company with local partners to produce low volume, made-to-order instruments that are no longer offered by the original, European manufacturers.

## PT. HYGEIA MEDICAL INDONESIA

Formed to be our distributor in Indonesia, the plan is to emulate the Hymed concept and offer in-house warranty and repair services, before moving to finishing instruments, ultimately setting up a domestic manufacturing facility for the ASEAN market.

## HYGEIA INSTRUMENTS (MALAYSIA) SDN BHD

3a-4-139 Bukit Jambul Complex, Jalan Rumbia  
Bayan Lepas 11900 Pulau Pinang, Malaysia

The company is the authorised representative, sole importer and chief distributor of HygeiaMedical branded surgical instruments and a variety of other brands in Malaysia. The registered office of the company is in Selangor and principal place of operations is in Penang, with a sales office in KL.

The company is holding the MDA license number MDA-1590-WDP121 and is also authorised to participate in public projects and tenders with CIDB B0102202008250001070882.



Through its connections in the industry and due to its location in the major regional medical manufacturing hub of Penang, the company also acts as purchasing agent and supplies components and devices to customers in Europe, the United States and the Middle East.

## HYMED MANUFACTURING SDN BHD

The company was established in 2022 by Hygeia and local partners in Penang to provide local warranty services and repairs of surgical instruments. It is also planned to offer training facilities and technology transfer as part of a shift towards local manufacturing and assembly. One step towards this goal is the establishment of workshops in strategic locations with other local partners.

## HYPHARM HEALTHCARE SDN BHD

In 2023 Hygeia Medical and Hygeia Instruments established HYPHARM in order to expand into the IV testing and pharmaceutical sectors with focus on tropical diseases, common ailments and oncology. The company focuses on importing and distributing test kits and reputable generic medicines, but intends to expand into local and regional assembly and domestic production in the future.



### MDA Establishment License

EL-20210205-5242

Establishment Name : HYGEIA INSTRUMENTS (MALAYSIA) SDN BHD

License No : MDA-1590-WDP121

License Start Date : 12/01/2021

License End Date : 11/01/2024



### Good Distribution Practice for Medical Device - GDPMD

Establishment Name : HYGEIA INSTRUMENTS (MALAYSIA) SDN BHD

Certificate Number : MYG11206674

Original Issue Date : 16/11/2020

Recertification Date : 03/10/2023

Expiration Date : 15/11/2026



## MDA Registrations

GA1141121-60270 26/04/2023 25/04/2026	<b>SCISSORS</b>	Standard, operating, dressing, lister, supercut, micro, tungsten carbide, etc.	2020-04 2020-05 2020-06
GA1887621-61990 24/05/2021 23/05/2026	<b>DISSECTING INSTRUMENTS</b>	Scalpels, Autopay devices, Training instruments, etc.	2020-03
GA3246322-94832 17/06/2022 16/06/2027	<b>RETRACTORS</b>	All types, varieties and sizes Class A	2020-14
GA3334521-57808 24/03/2021 23/03/2026	<b>GENERAL MINOR SURGERY SET</b>	Selection of common instruments used in surgeries and clinics	MSSMy21
GA3779523-151033 19/10/2023 18/10/2028	<b>DISPOSABLE PROCTOSCOPE</b>	HM-PROCTOSCOPE	HM-PROC
GA4291721-76750 24/10/2021 23/10/2026	<b>UROLOGY</b>	Urological and Bladder instruments sets	2020-26
GA4298521-70964 10/09/2021 09/09/2026	<b>DIGESTIVE TRACT</b>	Intestines, Stomach, Rectum	2020-25
GA4318522-101689 22/08/2022 21/08/2027	<b>SUTURE &amp; DRESSING</b>	All types and sizes of needle holders, dressing forceps, etc.	2020-18
GA4677323-137989 30/06/2023 29/06/2028	<b>ARTERY COMPRESSES</b>	Class A - Artery Hemostasis Compression Devices	2023-AP ZXT/GF
GA7281923-133939 08/06/2023 07/06/2028	<b>LARYNGOSCOPE</b>	Blades, Handles, Sets & Packs Emergency & Surgical uses All sizes and shapes	2023-FO CO STD
GA7585023-128729 18/04/2023 17/04/2028	<b>DENTAL CHAIR</b>	Mikata brand Japanese technology electric powered chair for dental practice use	MKT/ U





**MDA Establishment License** EL-20231023-13137

Establishment Name : HYGEIA INSTRUMENTS (MALAYSIA) SDN BHD

License No : MDA-5254-WDP123

License Start Date : 12/01/2024

License End Date : 11/01/2027




EUDAMED SRN: IE-AR-000023941 Ireland Co. Reg. 618490




















**HYGEIA MEDICAL™ LIMITED**

Kyleteologue House, Tallanstown, Co.Louth, Ireland  
PF2513 Hindenburg Str., 46395 Bocholt, Germany

[www.hygeia.icu](http://www.hygeia.icu)




## MDA Registrations

GA8644422-92870 19/05/2022 18/05/2027	<b>CARDIO-VASCULAR</b>	Our complete Cardiovascular Surgery Set.	2020-21		
GA8861321-55543 22/02/2021 21/02/2026	<b>ESSENTIAL UTILITIES SET</b>	The most commonly used Scissors in the clinics and surgeries we supply	EUS1		
GA8647821-54377 04/02/2021 03/02/2026	<b>GENERAL HOSPITAL SET</b>	The most commonly used instruments in hospitals in Malaysia	GHSMY21		
GA10448021-65800 05/07/2021 04/07/2026	<b>GYNECOLOGY</b>	The extensive Gynecology catalogue is part of our OB/GYN section	2020-27		
GA10555421-65280 28/06/2021 27/06/2026	<b>OBSTETRICS</b>	The catalogue of Obstetrics devices is part of our OB/GYN section	2020-28		
GA10435223-130061 10/05/2023 09/05/2028	<b>OTOLOGY</b>	Our Otology catalogue is part of the ENT family of sets	2020-30		
GA10572423-134465 05/10/2023 04/10/2028	<b>RHINOLOGY</b>	Our Rhinology catalogue is part of the ENT family of sets	2020-31		
GA8363521-55129 17/02/2021 16/02/2026	<b>EYE CATARACT REMOVAL SET</b>	This set of fine micro-instruments is part of our ENT family of sets	ECRS21		
GA8143021-52273 21/02/2021 20/02/2026	<b>ORTHOPAEDIC SURGERY SET</b>	The most commonly used instruments in basic orthopedic surgeries.	MOSS21		
GA8819121-59520 15/04/2021 14/04/2026	<b>SCALPELS</b>	Blades, handles, specialties, all sizes and shapes	2020-01		
GA9297823-149650 05/10/2023 04/10/2028	<b>EXTERNAL FIXATOR SET</b>	We commissioned the large external fixator set after Aesculap discontinued their production.	LEFS23		
GC10606323-136059 21/06/2023 20/06/2028	<b>DISPOSABLE ENDOCUTTERS</b>	These disposable endocutters and their reusable handles are made for us by Lepu Medical	JQ/JJ		






## GENERAL HOSPITAL SET (MALAYSIA)

MDA Registration No.: GA8647821-54377

SCAN QR TO  
DOWNLOAD PDF








# BASIC ORTHOPAEDIC SURGERY SET (MALAYSIA)

MDA Registration No. GA8143021-55273

SCAN QR TO DOWNLOAD PDF








for soft wires  
Ø 2,8 mm  
for hard wires  
Ø 2,2 mm

**TCF-19-22**  
Wire Cutting  
TC Forceps  
(23cm)



**DF-198-2240**  
Standard Wire  
Pliers  
(16cm)



**DF-211-2563**  
McDonald Elevator  
double-ended  
(18.5cm)



**DF-197-2223**  
Boehler Extension  
Rod  
for femur and knee  
(210mm x 110mm)



**RD-34296-05**  
Orthopaedic Quick Coupling T Handle  
w. JACOBS 3-jaw-chuck, keyless locking  
opening up to 6.35mm Ø

**RD-34628-25**  
2.5mm 20cm

Orthopaedic Screwdriver  
with Ferrozell handle  
Hexagonal Head

**RD-34628-25**  
3.5mm 25cm



**STR-01-051**  
170mm  
1mm

**STR-01-052**  
170mm  
2mm

**STR-01-053**  
170mm  
3mm

**STR-01-054**  
170mm  
4mm

**Ceramic Bayonet Kerrison Ronguer**  
thin footplate, 40° FWD  
(easily disassembled for autoclave)



**GPC-913P-15**  
Stainless Steel Wire Rod  
plain/un-threaded 150mm

**GPC-913F-15**  
Stainless Steel Wire Rod  
fully threaded 150mm

**STR-01-026**  
150mm  
2mm

**STR-01-027**  
150mm  
3mm

**STR-01-028**  
150mm  
4mm

**STR-01-029**  
150mm  
5mm

**Footplate Kerrison Ronguer**  
150mm (thin) / 200mm (ultra thin)

**STR-01-030**  
200mm  
1mm

**STR-01-031**  
200mm  
2mm

**STR-01-032**  
200mm  
3mm

**Ceramic Bone Ejection Kerrison Ronguer**  
(Ceramic coated)  
200mm (thin footplate)



Stainless Steel Wire  
(316L at 10m length)  
**GPC-915-045** 0.45mm  
**GPC-915-055** 0.55mm  
**GPC-915-080** 0.80mm  
**GPC-915-100** 1.00mm  
**GPC-915-120** 1.20mm  
**GPC-915-150** 1.50mm

**STR-01-038**  
FWD Ejector  
2mm

**STR-01-039**  
FWD Ejector  
3mm

**STR-01-040**  
FWD Ejector  
4mm

**STR-01-041**  
FWD Ejector  
5mm

**STR-01-042**  
FWD Ejector  
6mm

**STR-01-043**  
No Ejector  
1mm

**STR-01-044**  
No Ejector  
1.5mm

\* stainless steel versions available under SS code, matte finish available with MS code

# CATARACT REMOVAL SET (GA8363521-55129)



SCAN QR TO  
DOWNLOAD PDF



**TCF-04-13**  
Castroviejo  
Needle  
Holder  
with lock  
(12.5cm)



**TCF-04-13**  
Barraquer  
Needle  
Holder  
(15cm)



**DF-371A-4331A**  
Castroviejo Scissors  
left (10cm)



**DF-371B-4315M**  
Castroviejo  
Corneal Scissors  
(10cm)


**DF-371B-4316P**  
Aebli Capsule Scissor:  
(10cm)



**DF-365-4266**  
Castroviejo Suture Forceps  
(10cm)



**DF-360-4214**  
Graeffe  
Forceps  
1x2 (7cm)



**DF-362-4242**  
Bishop-Harmann  
Fixation Forceps  
1x2 (8.5cm)

**DF-366-4278**  
Bonaccolto  
Suture Tying Platform  
smooth, angled  
(11cm)



**2018-EOL-M**  
Kelman-McPherson  
angled tying platform



**DF-4-132** Scalpel Handle No. 3 – solid  
(12cm)



**DF-362-4239**  
Hess Forceps  
(7cm)



**DF-362-4235**  
Gill-Hess Iris  
Forceps  
(9cm)

**DF-362-4234**  
Gill Iris Forceps  
(9cm)

**2019-EOL-M**  
Masket Capsulorheads  
very delicate sharp tips

**2020-EOL-M**  
Dardenne Suture Tying  
6mm tying platform  
(10cm)



**DF-341-389:**  
Iris knife  
(12cm)  
with 2 cutting  
edges  
of 5mm



**DF-348-402**  
Stevens  
Hook  
(12cm)



**DF-351-4074**  
Lewis  
Lens Scoop  
smooth  
(12cm)



**DF-350-4059**  
Barraquer  
Iris Spatula  
smooth (12cm)



**DF-350-4065**  
Castroviejo  
C Spatula  
(10cm)



**DF-354-4127**  
Wilder  
(10cm)




**DF-337-3818**  
Desmarres  
Lid Retractor

**DF-333-3778**  
Barraquer  
Adult size (4cm)


**DF-333-3777**  
Barraquer  
(Colibri)  
Child size (3cm)



**DF-350-4053**  
Castroviejo  
Spatula  
(12cm)



**DF-350-4049**  
Jaffe  
Lens Spatula  
(12cm)



**DF-350-4050**  
Kwitko Spatula  
(12cm)



**DF-350-4048**  
Hirschmann  
Spatula  
(12cm)



**FSI-SIC-23**  
Simcoe Irrigation Cannula  
23 gauges



**FSI-AIC-27**  
Air Injection Canula  
angled 5mm at blunt tip  
27 gauges



**FIS-RCC-21**  
Randolph Cykloidalysse Cannula  
angled 13mm from flat tip  
21 gauges



**FSI-LC-23**  
Lacrimal Cannula  
malleable tip  
23 gauges

Hygeia Medical also produces a large number of specialty instruments, sets and devices in cooperation with several well known brands from around the globe. Among these are dental implants, orthopedic implants, but also laparoscopic instruments for key-hole and minimal invasive procedures, many different types of alligators for ENT, as well as most kinds of laryngoscopes, electrosurgical devices, single use proctoscopes, devices and instruments for otology, ophthalmology and even personal blood glucose measuring devices.



We are currently offering two types of reusable (autoclavable) click-line-like handles with over 20 attachments, as well as a bipolar reusable forceps.

Also available are reusable clip applicators, needle holders, claw and spoon forceps, suture passers, trocars, retractors, grips, and suction components. Furthermore we have knot pushers, screws, needles, hooks and many more.

Lastly we offer disposable dissecting forceps, Metzenbaum scissors, and a variety of grasping forceps.

### Our Aesculap type **Large External Fixator Sets** (MDA Registration GA9297823-149650)

come either as fully compatible component sets, or as complete assembly kits, including all required tools and implementa for assembly and deployment. And we can get our workshops to tailor specific items to your particular need if and when required.


And the German-made orthopedic implants of the INTERCUS brand for whom Hygeia Instruments is the sole regional distributor are currently undergoing evaluation by the MDA.



We offer a number of Laryngoscopes for emergency and for surgical use in a variety of types, shapes & size ranges from neonatal all the way to XL adults. All blades and casings are made from top stainless steel and many of the blades use the latest in fiber optic technology, such as the Macintosh blades with adjustable tip and spare LED bulbs. All sets are packed in durable plastic containers with felt casing for safe storage and transportation.



## MEDICAL POWER TOOLS



TECH SPECS	S1000
Speed (rpm)	- 1,100
Torque (Nm)	> 3.0
Operation (V)	14.4V
Noise (db)	< 75

Our multi-functional medical electric drill sets can operate these attachment options:

- 1、Canulate drill.
- 2、Oscillating saw.
- 3、Bone drill
- 4、Dual function acetabulum reaming.drill
- 5、Cranial drill.
- 6、Cranial bur.
- 7、Sternum saw.
- 8、Reciprocating saw.
- 9、K-wire drill



## REPAIRS, WARRANTY SERVICE & MAINTENANCE HYMED MANUFACTURING SDN BHD

Through HyMed, a joint venture between Casselman Surgical in Penang and the Hygeia Group, we are able to offer in-country QC, warranty repairs and maintenance services for all products under the HygeiaMedical brand, typically in two weeks, or less. Furthermore HyMed can also repair and refurbish third party branded instruments and devices, quickly and efficiently

On this basis HyMed will also work closely with the upcoming Hygeia School of Medical Technology to establish training centres and workshops in strategic locations around the country.



If you are an importer or distributor who needs warranty services for foreign branded products, or a POC operator who needs instruments repaired and refurbished, or if you are looking for ways to better train your med techs and instrument maintenance staff, we are the people to talk to.

**We can help with contract manufacturing, regional component sourcing and local assembly.**

# HygeiaMedical



Hygeia Medical Ltd. Reg. 616490 EUDAMED SRN IE-AR-000023941 TRN 3957302JH

[www.hygeia.icu](http://www.hygeia.icu)



IRELAND - GERMANY - MALAYSIA - INDONESIA - PHILPPINES - AUSTRALIA